




Impact of high power and high incidence angles on peripheral prism for homonymous hemianopia




Jae-Hyun Jung and Eli Peli

Schepens Eye Research Institute, Massachusetts Eye and Ear, Harvard Medical School, Boston, MA



Introduction

- Homonymous Hemianopia (HH)** occurs due to stroke, head injury or brain surgery
- Peripheral Prisms**^{1,2} (2000)
Expand upper and lower segments of the lateral visual field using Fresnel prisms




Peli Prism glasses
Rigid Fresnel Peripheral Prisms

Binocular visual field of a person with HH

Binocular visual field expansion with 40Δ peripheral prisms

Prism Power Variation in High & Low Power Prisms



Distortion and Reflections in High Power (57Δ) Fresnel Prism

- Prism distortion**^{3,4}
Prism power variation results in **reduced apical scotoma** and **minification** (increased field expansion)
- Spurious reflections**
Dimly show unrelated regions of the scene (grid lines and blue blur from a window).
- Total Internal Reflection (TIR)**
Beyond the red arrow, only reflection portions (grid lines) are seen.

Prism Power Variation by Angle of Incidence

- Deflection angle is highly dependent on the **angle of incidence**
- Above a critical angle of incidence, **total internal reflection (TIR)**
- The size of **Apical scotoma** is same as the effective prism power at the apex.

OPS and EPS Configurations

- Eyeward Prism Serration (EPS):**
Configuration Commonly used for press-on prism
- Outward Prism Serration (OPS):**
Configuration commonly used for permanent PMMA Fresnel prism

Magnification and Minification in OPS and EPS

Simulated Views in Primary Gaze, 10° and 20° Gaze Shifts to the Blind Side

- With gaze shift to blind side**
EPS: Magnification is almost constant (almost CDA), minimal distortion
OPS: Increasing minification (distortion) until TIR, larger expansion dim and narrow

- At primary gaze**
EPS: Magnification, field expansion is smaller, apical scotoma is larger
OPS: Minification, field expansion is larger, apical scotoma is smaller

Spurious Reflections in 57Δ Fresnel prism

- Monocular diplopia, monocular visual confusion, false alarms, & contrast reduction** can be caused by spurious reflections
- Meaningful reflections (**forming images**) in Fresnel prism
- Light rays split on surfaces into **transmitted** and **surface reflected** rays

I: Intended image
TIR: TIR on base
SR: Surface reflection

Conclusion

- Prism power varies with angle of incidence
 - For low-power prisms (20Δ) the variations within the range of practical gaze shift (±15°) are very small
 - For high-power prisms (57Δ) The effective deflection varies with angle of incidence, and these effects must be considered when prescribing prisms for visual field expansion
- In **OPS** configuration
 - The field-of-view through the prism is **wider** and therefore more **compressed**
 - For 57Δ **TIR** blocks views beyond ≈5° into the blind side which **limits the effects of scanning**
 - Apical Scotoma** is smaller (can cause diplopia)
 - At primary gaze the **TIR on base** appears in the **blind hemifield** while **dimmer surface reflections** appear at the **seeing hemifield**
- In **EPS** configuration prism power and its variability are reduced
 - The actual prism view is slightly **magnified** by the reduced prism power. **The effect of scanning is not limited**
 - At primary gaze the **TIR on base** and the **surface reflection** appear in the **visible hemifield** and can cause **visual confusion, diplopia, & bright false alarms**

Acknowledgement and References

- This work was supported in part by NIH grants EY12890 and EY023385 (EP), and the Basic Science Research Program (2012R1A6A3A03038820), National Research Foundation of Korea (NRF), the Ministry of Education, Science and Technology (J-HJ)
 - Dr. Peli has patent rights (assigned to Schepens Eye Research Institute) for the peripheral oblique prisms (licensed to Chadwick Optical)
- E. Peli (2000)
Field expansion for homonymous hemianopia by optically-induced peripheral extropia, *Optometry and Vision Science*, **77**(9), 453-464
 - E. Peli (2008)
Peripheral Field Expansion Device, *United States patent* 7,374,284
 - H. L. Apfelbaum, N. C. Ross, A. B. Bowers, E. Peli (2013)
Considering optical scotomas, confusion, and diplopia when prescribing prisms for homonymous hemianopia, *Translational Vision Science & Technology* **2**(4), article 2
 - J. -H. Jung, E. Peli (2014)
Impact of high power and angle of incidence on prism corrections for visual field loss, *Optical Engineering* **53**(6), 061707, Open Access